

Validation Report

AlerTox ELISA Mustard

KIT3058/KT-5915

1.	Scope	2
2.	Precision	. 2
	A) Intra-Assay Variation	. 2
	B) Inter-Assay Variation	. 3
3.	Recovery	3
4.	Analytical Sensitivity	5
5.	Linearity	6
6.	Cross-Reactivity	. 8
7.	Robustness	. 9



1. Scope

The AlerTox ELISA Mustard is designed for the determination of Mustard in food. The present report describes the validation process and its results.

2. Precision

A) Intra-Assay Variation

The intra-assay variation was determined by testing three controls of various concentration levels in 20fold replicates.

Replicate	Level 1	Level 2	Level 3	
1	3.4	11.5	28	
2	3.8	12.8	32	
3	4.1	11.8	31	
4	4.7	12.2	32	
5	4.6	12.1	30	
6	4.3	12.2	31	
7	4.6	11.0	30	
8	4.6	13.5	30	
9	3.6	12.4	31	
10	3.8	12.5	29	
11	4.0	11.2	31	
12	4.3	12.3	31	
13	3.9	12.6	32	
14	3.8	12.3	29	
15	3.4	12.2	28	
16	4.2	13.2	33	
17	2.9	12.7	28	
18	3.5	12.2	32	
19	3.6	13.2	32	
20	5.1	13.6	35	
Mean	4.0	12.4	31	
SD	0.53	0.70	1.9	RMS
CV [%]	13.3	5.7	6.0	9.0

Table 1: Intra-assay variation of the AlerTox ELISA Mustard

The coefficient of variation is ranging from 5.7% to 13.3% depending on the concentration.



B) Inter-Assay Variation

The inter-assay variation was determined by testing three controls of various concentration levels in four different test runs of the same kit lot.

Assay No.	Level 1	Level 2	Level 3	
1	5.5	13.2	29	
2	4.7	17.2	30	
3	5.0	15.1	35	
4	4.4	20.0	34	
Mean	4.9	16.4	32	
SD	0.47	2.91	3.0	RMS
CV [%]	9.6	17.8	9.4	12.9

Table 2: Inter-assay variation of the AlerTox ELISA Mustard

The coefficient of variation is ranging from 9.4% to 17.8% depending on the concentration.

3. Recovery

For recovery experiments different sample matrices were spiked with mustard to obtain various final concentrations after performing all sample pre-treatment steps. Tested samples and results were as follows.

Table 3: Recovery of various samples tested with the AlerTox ELISA Mustard

Canned Soup					
	Actual	Recovery [%]			
Target Value	Concentration				
5 ppm	5.5	110			
15 ppm	13.4	89			
35 ppm	30.9	88			
	Mean	96			

Salad Sauce

Target Value	Actual Concentration	Recovery [%]
5 ppm	4.1	81
15 ppm	10.2	68
35 ppm	28.0	80
	Mean	76



Cheese

	Actual	Recovery [%]
Target Value	Concentration	
5 ppm	4.3	86
15 ppm	11.7	78
35 ppm	35.8	102
	Mean	89

Sausage

	Actual	Recovery [%]
Target Value	Concentration	
5 ppm	5.1	102
15 ppm	14.9	99
35 ppm	32.7	93
	Mean	98

Instant Soup

Target Value	Actual Concentration	Recovery [%]
5 ppm	5.0	101
15 ppm	10.6	71
35 ppm	24.2	69
	Mean	80

Mixed Herbs

	Actual	Recovery [%]
Target Value	Concentration	
5 ppm	3.9	77
15 ppm	11.3	75
35 ppm	28.7	82
	Mean	78

Mean recoveries are ranging from 76% to 98% depending on the sample matrix.



4. Analytical Sensitivity

For determination of the analytical sensitivity sample diluent was assayed in 24fold replicates. After identification of possible outliers the OD mean and standard deviation was calculated. The corresponding concentration of the OD mean + 3x standard deviation was defined as limit of detection. This results in limits of detection according to the following table:

Replicate	Sample diluent [OD]	
1	0.214	
2	0.212	
3	0.210	
4	0.198	
5	0.191	
6	0.185	
7	0.186	
8	0.185	
9	0.214	
10	0.189	
11	0.193	
12	0.170	
13	0.187	
14	0.176	
15	0.169	
16	0.157	
17	0.204	
18	0.209	
19	0.198	
20	0.179	
21	0.197	
22	0.181	
23	0.171	
24	0.166	
Mean	0.189	
SD	0.016	
Limit of Detection	1.3 ppm	

Table 4: Matrix-independent analytical sensitivity of the AlerTox ELISA Mustard

The limit of detection is 1.3 ppm of mustard. The lowest positive standard (2 ppm) was defined as limit of quantification (LOQ).



5. Linearity

Linearity was determined by spiking sausage, salad sauce, canned soup, instant soup, cheese and mixed herbs samples with mustard and testing subsequent dilutions of the resulting extracts. For calculation of the linearity the highest concentration was defined as reference value (100%) and further dilutions were expressed in percent of this reference after consideration of the dilution factor

Table 5: Matrix dependent linearity of the AlerTox ELISA Mustard

Canned Soup				
Target Value	Concentration [ppm]	Recovery [%]		
60 ppm	60.1	100		
30 ppm	31.7	105		
15 ppm	15.1	100		
7.5 ppm	7.2	96		
3.75 ppm	3.5	93		
	Mean [%]	99		

Salad Sauce

	Concentration	
Target Value	[ppm]	Recovery [%]
60 ppm	43.8	100
30 ppm	21.2	97
15 ppm	11.2	102
7.5 ppm	5.7	104
3.75 ppm	2.5	91
	Mean [%]	99

Sausage

	Concentration	
Target Value	[ppm]	Recovery [%]
60 ppm	55.6	100
30 ppm	26.3	95
15 ppm	13.0	94
7.5 ppm	6.3	91
3.75 ppm	3.1	89
	Mean [%]	92



Cheese

Target Value	Concentration [ppm]	Recovery [%]
60 ppm	55.7	100
30 ppm	27.7	99
15 ppm	12.2	88
7.5 ppm	6.9	99
3.75 ppm	2.8	80
	Mean [%]	92

Instant Soup				
Target Value	Concentration [ppm]	Recovery [%]		
60 ppm	43.7	100		
30 ppm	23.5	108		
15 ppm	12.8	117		
7.5 ppm	6.5	119		
3.75 ppm	2.9	106		
	Mean [%]	113		

Mixed Herbs

Target Value	Concentration [ppm]	Recovery [%]
60 ppm	53.0	100
30 ppm	26.7	101
15 ppm	14.2	107
7.5 ppm	6.6	100
3.75 ppm	3.3	100
	Mean [%]	102

For different matrices the mean linearity is ranging from 92% to 113%. The linearity is independent of the specific concentration and may only be affected by the intra-assay and inter-assay variation.



6. Cross-Reactivity

For the following foods no cross-reactivity (results < LOQ) could be detected:

Adzuki	Cow's milk	Peanut	
Almond	Cumin	Pecan	
Apricot	Curcuma	Pepper	
Barley	Dill	Pine seed	
Bean, white	Duck	Pistachio	
Beef	Fennel	Poppy seed	
Bovine gelatin	Fenugreek	Pork	
Brazil nut	Garden cress	Potato	
Buckwheat	Garlic	Pumpkin seed	
Caraway	Gliadin	Radish	
Cardamon	Goat's milk	Rice	
Carob gum	Guar gum	Rye	
Carrot	Hazelnut	Saccharose	
Cashew	Horseradish	Sesame	
Cayenne	Kidney bean	Shrimp, cooked	
Celery	Kiwi	Shrimp, raw	
Cherry	Lamb	Shrimps	
Chestnut	Leek	Soy flour	
Chia	Lentil	Soy lecithin	
Chicken	Linseed	Split peas	
Chickpea	Lupin	Sunflower seeds	
Chili	Macadamia	Thyme	
Cinnamon	Nutmeg	Tomato	
Clove	Oats	Turkey	
Сосоа	Onion	Walnut	
Coconut	Paprika	Wheat	
Cod	Реа	White cabbage	
Corn	Peach		

Table 6: Non-cross-reactive food matrices in the AlerTox ELISA Mustard

The following cross-reactivity could be determined:

Food	Cross-reactivity [%]
Brown mustard (seed)	59
Black mustard (seed)	50
Charlock mustard	48
Rapeseed	59



7. Robustness

Robustness was determined by variation of different handling parameters as defined in the instruction manual. The results were compared with the results of samples analyzed according to the intended method. An un-spiked soup sample and a sample spiked with 15 ppm of mustard were analyzed respectively.

Variation of extraction temperature

The extraction temperature, defined as 60 °C, was changed to 25 °C, 40 °C and 70 °C, respectively.

Table 8: Variation of extraction temperature in the AlerTox ELISA Mustard

Sample	Result 60 °C	Result 25 °C	Result 40 °C	Result 70 °C
Soup 0 ppm	0 ppm	0 ppm	0 ppm	0 ppm
Soup 15 ppm	11.9 ppm	12.0 ppm	11.4 ppm	11.2 ppm

Under consideration of the intra-assay and inter-assay variations as stated in chapter 2 the results do not differ significantly.

Variation of extraction time

The extraction time, defined as 15 min, was changed to 5 min. 10 min and 20 min. respectively.

Table 9: Variation of extraction time in the AlerTox ELISA Mustard

Sample	Result 15 min	Result 10 min	Result 20 min
Soup 0 ppm	0 ppm	0 ppm	0 ppm
Soup 15 ppm	11.9 ppm	11.0 ppm	10.7 ppm

Under consideration of the intra-assay and inter-assay variation as stated in chapter 2, the results do not differ significantly.

Drift

In contrast to the test procedure as defined in the instruction manual the incubation time of the samples was extended and reduced by 4 minutes compared to the calibrators (20 min).

Table 10: Drift in the AlerTox ELISA Mustard

Sample	Result 20 min	Result 16 min	Result 24 min
Soup 0 ppm	0 ppm	0 ppm	0.2 ppm
Soup 15 ppm	11.9 ppm	9.8 ppm	13.5 ppm

The results differ significantly. Drift in extensive test runs should be avoided by pipetting calibrators once before the samples and once after the samples, using the mean value for calculation.

Date: Jul 2018