Detection of Salmonella in Raw and Ready-to-Cook Meat and Poultry Products and Environmental Samples Using an ISO-Validated Multiplex PCR Kit: the foodproof[®] Salmonella Genus plus Enteritidis & Typhimurium Detection LyoKit

Identification of the most relevant serovars for human and animal health, Salmonella Enteritidis and Salmonella Typhimurium, is critical when managing the food production chain from farm-to-fork to reduce consumer risk. Global regulations concerning these serovars require test methods to be robust in sensitivity and specificity to ensure reliable and accurate results.

Validation Methods

- The study comprised 75 Salmonella spp., 25 Salmonella Enteritidis and 25 Salmonella Typhimurium inclusivity strains, 30 exclusivity non-target strains, sensitivity, relative level of detection (RLOD) and collaborative studies.
- Validation of the alternative method was performed and compared to the ISO reference method (culture) for: .
 - Raw Meat and Poultry (25 g) 0
 - Ready-to-Cook Meat and Poultry (25 g) 0
 - 0 Environmental Samples
- Samples were enriched in:
 - Buffered Peptone Water (1:10) at 37 ± 1 °C for 16 20 h 0
- Following incubation, DNA extraction was performed with • foodproof[®] StarPrep[®] Three Kit (Standard or 8-strip), then lysates were analyzed by real-time PCR with the foodproof Salmonella Genus plus Enteritidis & Typhimurium Detection LyoKit.

Validation Results

- The validation study (NordVal No. 055) indicated that the foodproof method performs comparably to reference methods EN ISO 6579-1:2017 and 6579-3:2014 and fulfills the validation criteria according to EN ISO 16140-2:2016.
- The specificity studies yielded:
 - 100% inclusivity of all target strains 0
 - 100% exclusivity of all non-target strains 0
- Additionally, the acceptability limits for the sensitivity and RLOD studies for all categories and enrichment protocols were met.

Industry Significance

- The validation of this foodproof multiplex, real-time PCR assay provides many industries with a rapid and reliable method for the detection of Salmonella spp., Salmonella Enteritidis and Salmonella Typhimurium in raw and ready-to-cook meat and poultry products and environmental samples.
- The foodproof assays have a wide range of extraction options and compatible instruments, including the • BAX[®] Q7 System and the Roche LightCycler[®] 480 II.
- Using a single assay to not only screen for Salmonella but also identify regulated strains with the same enrichment, lysate, and assay reduces costs and improves operational efficiencies throughout food production and environmental monitoring, as well as for high-throughput testing laboratories.





High-Throughput 8-Strip

Extraction Option









Standard Extraction Option





Products

Product No.	Description	Quantity
Enrichment Options		
MED2010	Buffered Peptone Water	2.5 kg
MED2011		500 g
DNA Extraction Options		
KIT230187	foodproof [®] StarPrep [®] Three Kit, Standard	96 reactions
KIT230188	foodproof [®] StarPrep [®] Three Kit, 8-Strip	480 reactions
PCR Test Options		
KIT230134 (LP) KIT230135 (RP) KIT230136 (DP)	foodproof [®] <i>Salmonella</i> Genus plus Enteritidis & Typhimurium Detection LyoKit	96 reactions
Real-Time PCR Instrument Options		
Multiple packages*	BAX [®] Q7 System	1 system

* Contact us at <u>www.hygiena.com/contact</u>

BAX[®], foodproof[®], Hygiena[®] and StarPrep[®] are registered trademarks of Hygiena. LightCycler[®] is a registered trademark of Roche Diagnostics.