

AlerTox® Sticks Walnut

Validation Report

Cat No.	KIT3092 (10 tests)
Test Name	AlerTox Sticks Walnut
Description	Immunochromatographic rapid test for qualitative determination of walnut antigen in foods, kitchen and production facilities.
Specificity and Sensitivity	<p>The LOD (limit of detection) of AlerTox Sticks Walnut is 2.25 ppm of walnut protein. The range of detection (ROD) is 2.25-100,000 ppm of walnut protein. Above this range, the test can present a hook effect.</p> <p>AlerTox Sticks Walnut uses a combination of monoclonal antibodies against a major English/common walnut antigen, the 11 S globulin seed storage protein, known as allergen Jug r 4 from <i>Juglans regia</i>. AlerTox Sticks Walnut does NOT detect the antigens of cereals, legumes and other nuts, including pecan, black walnut, peanut, hazelnut, almond, macadamia, pistachio, cashew, brazil nut and coconut.</p>
Storage	Store at 10-30 °C (50-86 °F)

Assay Procedure

Solid, liquid and surface samples: according to manual of use INS3027

Validation Protocol

35 samples were analyzed separately to determine the correlation between AlerTox Sticks Walnut and ELISA against Walnut antigen. The following data was obtained:

	Sample	Brand	ELISA Test	AlerTox Sticks Walnut
1	Protein bar with Hazelnut	012 Nutrition	negative	negative
2	Protein bar with Almond	012 Nutrition	negative	negative
3	Protein cookie coconut	Royal cake	positive	positive
4	Protein bar with almond	Bombbar	negative	negative
5	Chocolate bar with walnut	Alpen Gold	positive	positive
6	Milk chocolate with whole almonds	Milka	negative	negative
7	Mars	Mars, Inc.	negative	negative
8	Chocolate bar	Alpen Gold	negative	negative
9	Geisha chocolate with hazelnut filling	Fazer	negative	negative
10	Milk Chocolate	Karl Fazer	negative	negative
11	Milk chocolate with hazelnuts	Dove	negative	negative
12	Hazelnut drink	ALPRO	negative	negative
13	Walnut drink	OraSi	positive	positive
14	Soy drink	ALPRO	negative	negative
15	Soy calcium drink	NATRUE	negative	negative
16	Pistachio drink	Simple Foods	negative	negative
17	Cashew drink	ALPRO	negative	negative
18	Pistachio almond ice cream	Baskin Robbins	negative	negative
19	Cheese with nuts	Cheese Lovers	positive	positive
20	Granola bar with nuts and honey	Everyday	positive	positive
21	Swiss roll	7 days	negative	negative
22	Soft cakes cherry	Fine Life	negative	negative

	Sample	Brand	ELISA Test	AlerTox Sticks Walnut
23	Biscuits filled with lemon cream	Grisbi	negative	negative
24	Napolitaner original	Manner	negative	negative
25	Dark chocolate and hazelnut cookies	Merba	negative	negative
26	Cashew coconut butter	Natty's	negative	negative
27	Marzipan	ZENTIS	negative	negative
28	Nutella	Ferrero	negative	negative
29	Italian bread with black sesame, bran and olive oil	Baker House	negative	negative
30	Milkshake with the taste of peanut ice cream and nougat	DANONE	negative	negative
31	Yogurt with nuts, chia seeds and vanilla	DANONE	positive	positive
32	Dairy dessert with walnuts	Hacendado	positive	positive
33	Raw walnuts in rice flour	Hygiena reference material	positive	positive
34	Raw walnut in wheat flour	Hygiena reference material	positive	positive
35	Raw black walnut	nut.com	positive (low)	negative

According to these results, the correlation between both techniques is 100% in the samples analyzed.

Surfaces:

A surface of stainless steel was thoroughly cleaned with water and soap, then rinsed with 60% ethanol, washed again with water and soap and finally rinsed with distilled water and left to dry.

Test areas of 4x4 cm (16cm²) were drawn on the surface and spiked with different levels of walnut extract. The amount of sample was dispensed with a micropipette and distributed throughout the whole surface with the micropipette tip. The areas were left to dry at room temperature, covered to avoid dust.

Spiking samples are prepared from a raw walnut extract and dilutions thereof. Spikings were made at the concentrations shown in the following table (assuming a protein content of 15 g/100 g). Each concentration was tested in duplicate.

Amount of walnut in extract (ppm)	Amount of walnut in test area ($\mu\text{g}/16 \text{ cm}^2$)	Approximate amount of walnut protein in test area ($\mu\text{g}/16 \text{ cm}^2$)	Result
2500	100	15	Positive
2500	70	10.5	Positive
2500	50	7.5	Positive
1250	25	3.75	Positive
625	12.5	1.9	Positive
312	10	1.5	Negative
312	6.25	0.9	Negative
156	3.12	0.5	Negative
78	2.3	0.34	Negative
39	1.56	0.23	Negative

A faint positive was detected at $1.9 \mu\text{g}/16 \text{ cm}^2$ of walnut protein and therefore the limit of detection of the swab protocol was set to $2 \mu\text{g}/16 \text{ cm}^2$.