# 5-in-1 Water Test

REF KIT3063 (KT-6372) / KIT3064 (KT-6123)

INS3023 REV A

## 1. Test principle

Enzymatic test strips are intended for quick visual semi-quantitative determination of total hardness, pH, nitrites, nitrates and free chlorine in water. One test strip is intended for testing of a single water sample. **Determination of total water hardness** is based on a chemical reaction between an indicator dye – eriochrome black – with  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  cations, which causes concentration-dependent change of color of the indicator dye.

**Determination of pH** is based on a color reaction of a mixed acidic-alkaline indicator, the color is dependent on concentration of  $\text{H}^+$  cations. **For determination of nitrites ( $\text{NO}_2^-$ )**, a modified Griess reaction is used; the chromogen changes its color intensity depending on nitrites concentration.

**Determination of nitrates ( $\text{NO}_3^-$ )** is based on enzymatic reduction of  $\text{NO}_3^-$  to  $\text{NO}_2^-$  followed by the above mentioned modified Griess reaction. **Free chlorine determination** is based on TMB oxidation by chlorine leading to formation of a colored product with color intensity depending on chlorine concentration in water. Concentration of each component in water is estimated by comparison of color intensity of the corresponding sensor pad on a strip with its counterpart on the color scale.

## 2. Test sensitivity

- Total hardness (dH). Determination range:  $4^\circ\text{d} - \geq 16^\circ\text{d}$  ( $1^\circ\text{d} = 0.3566 \text{ mg} - \text{eq/l}$ ). Sensitivity:  $4^\circ\text{d}$  (1.5 mg-eq/l).
- pH. Determination range: pH 5 -  $\geq 9$ . Accuracy: 0.5 pH units.
- Nitrites ( $\text{NO}_2^-$ ). Determination range for nitrites: 1 mg/l -  $\geq 10$  mg/l. Sensitivity: 1 mg/l.
- Nitrates ( $\text{NO}_3^-$ ). Determination range for nitrates: 10 mg/l -  $\geq 250$  mg/l. Sensitivity: 10 mg/l.
- Free chlorine (Cl) Determination range: 1 mg/l -  $\geq 20$  mg/l. Sensitivity: 1 mg/l.

## 3. Contents

- Plastic vial with 25 or 50 test strips.
- Instructions for use.

## 4. Samples

Water from pools, aquariums, wells, tap water and any other open or closed water sources.

## 5. Sample preparation

**Before testing, all samples should be brought to room temperature 15 - 25 °C (59 - 77 °F); low sample temperature causes a decrease in test sensitivity; testing of hot samples is not possible!** Water sample should be taken into a clean container. Testing in the flow is not possible.

## 6. Test procedure

A simple timer and a well-lit place are required to perform the test.

1. All test components and samples should be brought to room temperature 15 - 25 °C (59 - 77 °F) before running the test.


2. Take the required number of test strips from the vial and close it IMMEDIATELY to avoid exposure to humidity.

3. Dip the strip into the water sample. All the sensor pads should be submerged.

4. After 2-3 seconds, take the strip out and remove excess of water by shaking it off or by careful touching the strip edge with a clean filter paper (paper tissue).

5. Place the test strip on an even, clean, dry surface with sensor pads up. After two minutes compare the color of each sensor pad with its corresponding zone on the color scale\*:

Normal ranges for each parameter are determined by state or local sanitary authorities.

\*The range limits shown on the color scale are for reference only.

## Important notes

- The test strips should be stored in a dry place at 15 - 25 °C (59 - 77 °F). Keep away from humidity, vapors of acids, alkali and organic solvents.
- Avoid direct sunlight exposure of the test strips and the color scale.
- Color intensity of color scales from different lots may vary – please refer only to the scale supplied with the same lot.
- Test strips should be used NO LATER THAN 15 minutes after removing them from the vial.
- After the first opening of the vial, all test strips should be used within 6 months.
- All the components of the test kit are disposable; do not re-use them.
- DO NOT TOUCH the sensor pads.
- Do not use the strips if the tube is obviously broken and/or damaged.
- Do not use the test strips beyond the expiration date.

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