



## AlerTox<sup>®</sup> ELISA Almond Kit

### Claim Support Summary

**REF** KIT3049





## Table of Contents

1. Scope.....	3
2. Precision.....	3
2.1 Intra-Assay Variation .....	3
2.2 Inter-Assay Variation .....	6
2.3 Reproducibility of Extraction .....	6
3. Analytical Sensitivity .....	7
4. Recovery.....	9
4.1 Food Samples.....	9
4.2 Heat-treated Matrices .....	12
4.3 Rinse and Clean-in-Place (CIP) Water .....	12
4.4 Swabs .....	14
5. Linearity .....	14
6. Cross-Reactivity.....	17
7. Biological or Process Variants .....	18
8. Robustness.....	18
8.1 Variation of Extraction Temperature.....	18
8.2 Variation of Extraction Time .....	19
8.3 Drift.....	19
9. Summary .....	19



## 1. Scope

The AlerTox® ELISA Almond Kit is based on a polyclonal antibody and designed to determine almond content in food. This report describes the validation process and its results.

### Definitions

CV = Coefficient of variation

IQR = Interquartile range

LOD = Limit of detection

LOQ = Limit of quantification

OD = Optical density

ppm = mg of almond protein per kg or L of sample (mg/kg or mg/L)

PT = Proficiency testing

SD = Standard deviation

## 2. Precision

### 2.1 Intra-Assay Variation

The intra-assay variation (Repeatability) was determined by testing a control at 3 concentrations (0.4, 1 and 8 ppm) in 24-fold replicates, using either sample diluent (Table 1) or a cookie matrix (Table 2). After identification of possible outliers, the following results were obtained:

**Table 1. AlerTox ELISA Almond Kit Intra-Assay Variation with Sample Diluent.**

Replicate	0.4 ppm	1 ppm	8 ppm
1	0.32	1.00	8.27
2	0.33	1.06	8.80
3	0.36	0.96	8.84
4	0.35	1.04	8.67
5	0.37	0.99	8.60
6	0.37	0.98	8.35
7	0.39	0.96	7.38
8	0.37	0.99	8.11
9	0.37	0.99	8.72
10	0.43	0.97	8.12
11	0.39	0.97	7.97
12	0.41	0.98	8.44
13	0.44	0.97	7.39
14	0.44	0.98	8.08
15	0.37	1.01	8.03
16	0.39	1.02	8.53
17	0.43	0.96	7.89
18	0.44	1.06	7.33
19	0.44	1.07	7.52
20	0.41	1.07	7.93
21	0.32	1.04	7.70
22	0.37	1.06	8.64
23	0.41	0.99	8.04
24	0.39	1.00	7.02
Mean	0.39	1.00	8.10
SD	0.04	0.04	0.51
CV [%]	9.4	3.8	6.3

Mean
6.5

The coefficient of variation (CV) with sample diluent ranged from 3.8 to 9.4%, depending on the concentration tested.

**Table 2. AlerTox ELISA Almond Kit Intra-Assay Variation with a Cookie Matrix.**

Replicate	0.4 ppm	1 ppm	8 ppm	
1	0.43	0.96	8.2	
2	0.49	0.92	8.3	
3	0.47	0.99	8.3	
4	0.56	0.92	7.7	
5	0.52	0.99	9.0	
6	0.54	0.88	8.9	
7	0.49	0.92	8.8	
8	0.45	0.96	8.3	
9	0.44	0.97	7.7	
10	0.50	0.95	8.5	
11	0.47	0.99	8.0	
12	0.56	0.93	7.9	
13	0.51	0.94	8.8	
14	0.47	0.95	8.6	
15	0.44	0.92	8.2	
16	0.46	0.95	7.9	
17	0.47	0.98	7.9	
18	0.55	1.00	8.5	
19	0.53	0.95	8.9	
20	0.51	0.92	9.0	
21	0.46	0.86	8.7	
22	0.47	0.86	8.3	
23	0.49	0.88	8.1	
24	0.52	0.95	8.7	
Mean	0.49	0.94	8.4	
SD	0.04	0.04	0.41	Mean
CV [%]	8.0	4.1	4.9	5.7

The CV with a cookie matrix ranged from 4.1 to 8.0%, depending on the concentration tested.



## 2.2 Inter-Assay Variation

The inter-assay variation (Reproducibility) was determined by testing a control at 3 concentrations (0.4, 1 and 8 ppm) in 4 different test runs of the same kit lot, using either sample diluent (Table 3) or a cookie matrix (Table 4).

**Table 3. AlerTox ELISA Almond Kit Inter-Assay Variation with Sample Diluent.**

Assay Number	0.4 ppm	1 ppm	8 ppm	
1	0.38	0.99	8.98	
2	0.43	1.05	8.64	
3	0.48	0.92	8.10	
4	0.49	0.90	8.34	
Mean	0.44	0.96	8.51	
SD	0.05	0.07	0.38	Mean
CV [%]	<b>10.6</b>	<b>6.9</b>	<b>4.5</b>	<b>7.3</b>

The CV with sample diluent ranged from 4.5 to 10.6%, depending on the concentration tested.

**Table 4. AlerTox ELISA Almond Kit Inter-Assay Variation with a Cookie Matrix.**

Assay Number	0.4 ppm	1 ppm	8 ppm	
1	0.40	0.95	8.1	
2	0.38	1.01	8.1	
3	0.35	0.99	8.8	
4	0.41	0.97	8.8	
Mean	0.39	0.98	8.5	
SD	0.03	0.02	0.41	Mean
CV [%]	<b>6.7</b>	<b>2.4</b>	<b>4.8</b>	<b>4.7</b>

The CV with a cookie matrix ranged from 2.4 to 6.7%, depending on the concentration tested.

## 2.3 Reproducibility of Extraction

To determine the reproducibility of the extraction process, an almond proficiency test (PT) material (DLA 14/2019, sample 1) was analyzed. The PT material provider assessed the material's homogeneity by analyzing 20 different samples with 3 different test kits. According to the PT report, the micro tracer analysis of the material resulted in a probability of 85%.

The material was extracted 10 times, and the CV was calculated using the mean results of all extractions (Table 5).

**Table 5. Reproducibility of Extraction of the AlerTox ELISA Almond Kit.**

Extraction Number	Mean Result [ppm]
1	33.2
2	37.1
3	32.4
4	30.4
5	38.5
6	32.7
7	33.2
8	34.1
9	32.2
10	35.1
Mean	33.9
SD	2.42
CV [%]	<b>7.2</b>

The CV was 7.2 for the inter-extraction reproducibility results.

### 3. Analytical Sensitivity

To determine the analytical sensitivity of the AlerTox ELISA Almond Kit, sample diluent and various representative, almond-free sample matrices were assayed in 24-fold replicates. After identification of possible outliers, the mean and standard deviation of the optical density (OD) were calculated. The corresponding concentration of the mean OD + 3X standard deviation (SD) was defined as the limit of detection (LOD) (Tables 6a and 6b).

**Note:** In some cases, various samples of the same kind (e.g., different cookies) were tested. Table 6a represents the matrix with the highest LOD, while Table 6b summarizes the results.

**Table 6a. Matrix-Independent and Matrix-Dependent Analytical Sensitivity of the AlerTox ELISA Almond Kit (Complete Data).**

Replicate	Sample Diluent [OD]	Cookie Matrix [OD]	Sweets Matrix [OD]	Ice Cream Matrix [OD]	Chocolate Matrix* [OD]	Muesli Matrix [OD]	Food Supplement Matrix [OD]
1	0.033	0.061	0.075	0.059	0.128	0.043	0.025
2	0.036	0.064	0.082	0.066	0.128	0.042	0.026
3	0.038	0.068	0.068	0.063	0.128	0.045	0.024
4	0.040	0.051	0.085	0.051	0.120	0.041	0.021
5	0.041	0.059	0.066	0.066	0.118	0.039	0.023
6	0.046	0.063	0.058	0.051	0.130	0.042	0.028
7	0.034	0.060	0.075	0.050	0.123	0.049	0.028
8	0.048	0.053	0.063	0.048	0.119	0.035	0.033
9	0.040	0.054	0.060	0.053	0.124	0.043	0.022
10	0.034	0.051	0.075	0.043	0.124	0.046	0.024
11	0.034	0.055	0.067	0.045	0.112	0.041	0.019
12	0.047	0.058	0.075	0.056	0.131	0.035	0.020
13	0.036	0.060	0.059	0.043	0.115	0.051 <sup>†</sup>	0.019
14	0.041	0.060	0.068	0.048	0.125	0.042	0.023
15	0.038	0.057	0.073	0.042	0.123	0.040	0.031
16	0.045	0.056	0.080	0.046	0.135	0.036	0.030
17	0.037	0.058	0.059	0.055	0.128	0.039	0.020
18	0.040	0.062	0.083	0.045	0.120	0.036	0.017
19	0.043	0.051	0.062	0.048	0.116	0.050	0.019
20	0.043	0.052	0.085	0.069	0.130	0.042	0.022
21	0.046	0.048	0.078	0.044	0.115	0.040	0.019
22	0.046	0.048	0.084	0.057	0.117	0.035	0.022
23	0.033	0.048	0.070	0.055	0.121	0.036	0.026
24	0.037	0.060	0.079	0.047	0.129	0.041	0.025
Mean	0.040	0.057	0.072	0.052	0.123	0.041	0.024
SD	0.005	0.005	0.009	0.008	0.006	0.005	0.004
LOD [ppm]	<b>0.06</b>	<b>0.07</b>	<b>0.06</b>	<b>0.02</b>	<b>0.20</b>	<b>&lt; 0.005</b>	<b>&lt; 0.005</b>
Corresponding Calibration Curve Data [OD]							
<b>0 ppm</b>	0.040	0.046	0.069	0.068	0.055	0.052	0.039
<b>0.4 ppm</b>	0.144	0.195	0.260	0.233	0.229	0.262	0.170
<b>1 ppm</b>	0.332	0.350	0.532	0.459	0.467	0.414	0.363
<b>4 ppm</b>	1.037	1.148	1.482	1.364	1.498	1.410	1.038
<b>10 ppm</b>	1.910	2.005	2.692	2.387	2.455	2.536	1.934

\* Used the AlerTox Polyphenol Additive during sample extraction.

† Identified as an outlier by 1.5 interquartile range (IQR) testing and not considered in the calculations.

**Table 6b. Matrix-Independent and Matrix-Dependent Analytical Sensitivity of the AlerTox ELISA Almond Kit (Summary).**

Matrix Type	Variant	LOD (ppm)
Sample diluent	N/A	0.06
Cookies	Butter cookie	0.07
	Oat-spelt cookie	0.05
	TUC cookie	0.03
Sweets	Type 1	0.01
	Candy bar	0.06
	Vitamin drop	0.01
Ice cream	Type 1	0.08
	Type 2	< 0.005
Chocolate*	Yogurt cubes	0.05
	Sarotti, 72% cocoa	0.20
	Lindt, 70% cocoa	0.05
Muesli	Type 1	< 0.005
Food supplement	Lab sample	0.01
	Almased	< 0.005
	Wellmix Classic	< 0.005

\* Used the AlerTox Polyphenol Additive during sample extraction.

LODs vary from < 0.005 to 0.20 ppm, depending on the sample matrix. Note that the derived LODs are strictly dependent on the CV and thus may vary in individual tests.

The lowest positive standard was defined as the limit of quantification (LOQ) to ensure that measurements of all uncontaminated matrices result in concentrations lower than this value.

## 4. Recovery

### 4.1 Food Samples

For recovery experiments, different sample matrices were spiked with an almond extract at 3 final concentrations (0.4, 1 and 8 ppm) after performing all sample pre-treatment steps. The result of the unspiked matrix was subtracted before evaluation. The tested samples and calculated results are shown in Table 7.

**Table 7. Recovery of Various Samples Tested with the AlerTox ELISA Almond Kit.**

Cookie: Butter Cookie		Unspiked Matrix: 0.05 ppm
Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.45	113
1 ppm	1.02	102
8 ppm	8.11	101
	Mean	106

**Cookie: Oat-Spelt Cookie**

Unspiked Matrix: 0.04 ppm

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	1.13	113
1 ppm	4.16	104
8 ppm	19.7	98
Mean		105

**Cookie: TUC**

Unspiked Matrix: 0.03 ppm

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.37	94
1 ppm	0.92	92
8 ppm	8.08	101
Mean		96

**Muesli**

Unspiked Matrix: 0.00 ppm

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.48	119
1 ppm	1.24	124
8 ppm	8.31	104
Mean		115

**Sweets: Candy Bar**

Unspiked Matrix: 0.00 ppm

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.43	108
1 ppm	0.93	93
8 ppm	7.92	99
Mean		100

**Sweets: Vitamin Drop\***

Unspiked Matrix: 0.00 ppm

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.38	95
1 ppm	0.95	95
8 ppm	9.19	115
Mean		102

**Chocolate: Yogurt Cubes\***

Unspiked Matrix: 0.00 ppm

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.32	80
1 ppm	1.02	102
8 ppm	6.58	82
Mean		88

\* Used the AlerTox Polyphenol Additive during sample extraction.

**Chocolate: Sarotti (72% Cocoa)\*****Unspiked Matrix: 0.19 ppm**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.30	76
1 ppm	1.34	134
8 ppm	8.46	106
Mean		105

**Chocolate: Lindt (70% Cocoa)\*****Unspiked Matrix: 0.00 ppm**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.38	96
1 ppm	0.95	95
8 ppm	7.73	97
Mean		96

**Chocolate: Lab Sample\*****Unspiked Matrix: 0.00 ppm**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.29	71
1 ppm	0.72	72
8 ppm	5.68	71
Mean		71

**Food Supplement: Almased****Unspiked Matrix: 0.00 ppm**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.47	118
1 ppm	1.04	104
8 ppm	8.31	104
Mean		109

**Food Supplement: Wellmix Classic****Unspiked Matrix: 0.00 ppm**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.27	67
1 ppm	0.73	73
8 ppm	10.21	128
Mean		89

**Food Supplement: Wellmix Vanilla****Unspiked Matrix: 0.00 ppm**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.32	79
1 ppm	0.70	70
8 ppm	5.43	68
Mean		72

\* Used the AlerTox Polyphenol Additive during sample extraction.



Ice Cream: Yogurt		
Unspiked Matrix: 0.00 ppm		
Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.25	63
1 ppm	0.73	73
8 ppm	6.47	81
Mean		72

Ice Cream: Solero		
Unspiked Matrix: 0.00 ppm		
Spiked Value	Measured Concentration [ppm]	Recovery [%]
0.4 ppm	0.30	76
1 ppm	0.66	66
8 ppm	7.35	92
Mean		78

Mean recoveries ranged from 72 to 115%, depending on the sample matrix.

## 4.2 Heat-treated Matrices

To assess potential changes in the analytical sensitivity and recovery of processed food samples, baked or unbaked cookie matrices with or without almonds were analyzed. The spiked cookie matrix had 10 ppm of almond and was diluted by half down to 1.25 ppm. Each unspiked or spiked sample was then divided into 2 parts that were either left unbaked or baked at 150 °C for 20 min (Table 8). All tests were performed in duplicates.

**Table 8. Accuracy of the AlerTox ELISA Almond Kit with Heated-treated Matrices.**

Matrix	No Almond Spike [0 ppm]	Almond Spike [1.25 ppm]	Almond Spike [2.5 ppm]	Almond Spike [5 ppm]	Almond Spike [10 ppm]
Unbaked cookie	0.03	0.95	2.33	4.58	9.28
Baked cookie	0.00	0.92	2.21	4.48	9.07

No significant difference between the baked and unbaked cookie samples was observed.

## 4.3 Rinse and Clean-in-Place (CIP) Water

The AlerTox ELISA Almond Kit sensitivity with rinse and CIP water samples was determined by testing various concentrations of diluted almond extract in different commonly used, water-based rinse solutions. The pH was adjusted to  $8.2 \pm 0.5$  before running the test. Afterwards, the samples were diluted 1:5 in pre-diluted sample extraction buffer. The samples tested and calculated results are shown in Table 9.

**Table 9. Recovery with Rinse Solution Samples Tested Using the AlerTox ELISA Almond Kit.****Pure Water****Unspiked Matrix: 0.01 mg/L**

Spiked Value	Measured Concentration [mg/L]	Recovery [%]
0.1 mg/L	0.32	79
0.2 mg/L	0.76	95
0.4 mg/L	1.29	80
0.8 mg/L	2.76	86
1.6 mg/L	6.29	98
Mean		88

**0.5% Sodium Triphosphate****Unspiked Matrix: 0.06 mg/L**

Spiked Value	Measured Concentration [mg/L]	Recovery [%]
0.1 mg/L	0.35	87
0.2 mg/L	0.76	95
0.4 mg/L	1.27	79
0.8 mg/L	2.08	65
1.6 mg/L	5.50	86
Mean		82

**2% Sodium Hydroxide****Unspiked Matrix: 0.01 mg/L**

Spiked Value	Measured Concentration [mg/L]	Recovery [%]
0.1 mg/L	0.46	115
0.2 mg/L	1.01	126
0.4 mg/L	1.71	107
0.8 mg/L	2.78	87
1.6 mg/L	5.35	84
Mean		104

**0.5% Nitric Acid****Unspiked Matrix: 0.00 mg/L**

Spiked Value	Measured Concentration [mg/L]	Recovery [%]
0.1 mg/L	0.33	83
0.2 mg/L	0.88	109
0.4 mg/L	1.87	117
0.8 mg/L	2.84	89
1.6 mg/L	5.93	93
Mean		98



0.5% Citric Acid		Unspiked Matrix: 0.02 mg/L
Spiked Value	Measured Concentration [mg/L]	Recovery [%]
0.1 mg/L	0.46	114
0.2 mg/L	0.75	94
0.4 mg/L	1.89	118
0.8 mg/L	2.58	81
1.6 mg/L	5.65	88
Mean		99

Mean recoveries for rinse water ranged from 82 to 104%, depending on the kind of rinse solution.

#### 4.4 Swabs

To monitor swab performance, surfaces (5 cm x 5 cm) were contaminated with defined amounts of an almond-containing solution (0.8, 2.4 and 8 ng/cm<sup>2</sup>). Each amount was applied in duplicate, and the surfaces were allowed to dry completely.

The swabs were moistened in 1 mL pre-diluted sample extraction buffer. Then, the surfaces were swabbed thoroughly. The swab samples were redissolved in the initial 1 mL of buffer, and 100 µL of this solution was used in the ELISA. The mean results were calculated as the amount of almond per area (Table 10).

**Table 10. Recovery with Swabs Tested Using the AlerTox ELISA Almond Kit.**

Swabs		Unspiked Matrix: 0.00 ng/cm <sup>2</sup>	
Spiked Value	Measured Concentration [ng/cm <sup>2</sup> ]	CV [%]	Recovery [%]
0.8 ng/cm <sup>2</sup>	0.47	11.8	59
2.4 ng/cm <sup>2</sup>	1.24	6.9	52
8 ng/cm <sup>2</sup>	5.24	10.8	66
Mean		9.6	59

Recovery ranged from 52 to 66%. Mean recovery was 59%. Repeated swabbing of the same surface indicates that varying recoveries are an intrinsic feature of the method that can only be improved by multiple repetitions. Recoveries are affected by the solubility of proteins, residual liquid and complete swabbing of the surface area.

Note that surfaces were contaminated by protein extracts and not complete commodities. The latter are more particulate by nature and might be swabbed more easily. Thus, the approach chosen here reflects a worst-case scenario.

## 5. Linearity

Linearity was determined by spiking various matrices with an almond extract and testing serial dilutions (10 to 0.625 ppm) of the resulting matrices (Table 11). For calculating linearity, the highest concentration was defined as the reference value (100%), and further dilutions were expressed as a percent of this reference after accounting for the dilution factor. For the matrices extracted with the addition of the AlerTox Polyphenol Additive, further dilution was done with the standard extraction buffer that did not include the additive.



Table 11. Matrix-Dependent Linearity of the AlerTox ELISA Almond Kit.

**Cookie: Butter Cookie**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
10 ppm	10.3	Reference* (100)
5 ppm	4.8	100
2.5 ppm	2.44	96
1.25 ppm	1.29	90
0.625 ppm	0.62	101
	Mean	97

**Muesli**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
10 ppm	10.7	Reference* (100)
5 ppm	5.6	105
2.5 ppm	2.45	91
1.25 ppm	1.30	97
0.625 ppm	0.66	99
	Mean	98

**Chocolate: Lindt (70% Cocoa)<sup>†</sup>**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
10 ppm	9.6	Reference* (100)
5 ppm	4.8	101
2.5 ppm	2.46	102
1.25 ppm	1.18	98
0.625 ppm	0.60	100
	Mean	100

**Food Supplement: Almased**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
10 ppm	9.3	Reference* (100)
5 ppm	4.6	100
2.5 ppm	2.28	98
1.25 ppm	1.07	93
0.625 ppm	0.54	93
	Mean	96

\* The highest concentration was defined as the reference value (100%). Further dilutions were expressed as a percent of this reference after accounting for the dilution factor.

<sup>†</sup> Used the AlerTox Polyphenol Additive during sample extraction.

**Ice Cream: Solero**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
10 ppm	8.7	Reference* (100)
5 ppm	4.3	99
2.5 ppm	2.19	101
1.25 ppm	1.00	92
0.625 ppm	0.57	105
	Mean	99

**Sweets: Vitamin Drop<sup>†</sup>**

Spiked Value	Measured Concentration [ppm]	Recovery [%]
10 ppm	10.0	Reference* (100)
5 ppm	4.6	91
2.5 ppm	2.30	92
1.25 ppm	1.16	93
0.625 ppm	0.62	98
	Mean	93

\* The highest concentration was defined as the reference value (100%). Further dilutions were expressed as a percent of this reference after accounting for the dilution factor.

† Used the AlerTox Polyphenol Additive during sample extraction.

For different matrices, the mean linearity ranged from 93 to 100%. The linearity was independent of the specific almond concentration and may have been affected by the precision, as stated in Section 2.



## 6. Cross-Reactivity

No cross-reactivity is defined as a result that is less than the LOQ. For the AlerTox ELISA Almond Kit, the LOQ is 0.5 ppm whole almond. During testing, no cross-reactivity was detected in the foods included in Table 12.

**Table 12. Non-Cross-Reactive Food Matrices with the AlerTox ELISA Almond Kit.**

Non-Cross-Reactive Matrices			
Raw Material	Concentration [ppm]	Raw Material	Concentration [ppm]
Adzuki bean	0.11	Ginger (fresh)	0.00
Apricot	0.00	Gliadin	0.00
Bean, white	0.07	Guar gum	0.00
Beef	0.00	Gum arabic	0.02
Brazil nut	0.07	Hazelnut	0.10
Buckwheat	0.03	Horseradish	0.00
Cabbage, white	0.00	Milk, cow	0.04
Cardamom	0.05	Milk, goat	0.04
Carob gum	0.00	Pepper, black	0.12
Carrot	0.00	Pine seed	0.00
Cashew	0.00	Pistachio	0.03
Cayenne	0.05	Pork	0.01
Celery	0.27	Potato	0.06
Chestnut	0.28	Prawn	0.00
Chicken	0.00	Pumpkin seed	0.00
Chickpea	0.10	Radish	0.03
Chili	0.11	Rapeseed	0.31
Cinamon	0.09	Rice	0.00
Cocoa	0.39	Rye	0.23
Coconut	0.02	Saccharose	0.10
Cod	0.00	Sesame	0.04
Corn	0.00	Shrimps	0.00
Cumin	0.00	Soy flour	0.00
Dill	0.04	Soy lecithin	0.38
Duck	0.00	Split pea	0.00
Egg (dried)	0.05	Sunflower seed	0.00
Fennel	0.06	Thyme	0.03
Fenugreek	0.00	Tofu	0.00
Flaxseed	0.02	Tomato	0.00
Garden cress	0.00	Turkey	0.02
Garlic (fresh)	0.00	Turmeric	0.01
Garlic (granulated)	0.01	Walnut	0.14
Gelatin, cow	0.02	Wheat	0.22



The following cross-reactivities were determined (Table 13):

**Table 13. Cross-Reactive Food Matrices with the AlerTox ELISA Almond Kit.**

Cross-Reactive Food Matrices		
Raw Material	Concentration [ppm]	Cross-reactivity [%]
Apricot kernel	~1M	100
Peach kernel	~160,000	16
Cherry kernel	~17,000	1.7
Mahaleb cherry kernel	~14,000	1.4
Plum kernel	~10,000	1
Clove	1.0	0.0001
Poppy seed	0.7	0.00007
Chia	0.6	0.00006
Caraway	0.4	0.00004
Ginger (ground)	0.4	0.00004
Sweet red pepper seeds	0.4	0.00004

## 7. Biological or Process Variants

Industrial processing (e.g., heating, pH changes, fermentation) could result in structural variants of almond allergens that differ from the target allergen detected by the AlerTox ELISA Almond Kit. To assess the reactivity towards other variants, unroasted samples were compared to samples roasted at 150 °C for varying amounts of time (Table 14).

**Table 14. Accuracy (Relevant Analyte Variants) of the AlerTox ELISA Almond Kit.**

	Not Roasted	Roasted at 150 °C			
		15 min	30 min	60 min	120 min
Reactivity [%]	<b>100</b>	94	86	38	25

Variant-dependent reactivity ranged from 25 to 100%.

## 8. Robustness

Robustness was determined by the variation of different handling parameters defined in the instruction manual. The results obtained under various conditions were compared to the results obtained by following the instruction manual. An unspiked cookie sample and a sample spiked with almond extract were analyzed.

### 8.1 Variation of Extraction Temperature

The extraction temperature, defined as 60 °C, was changed to 25, 40 and 70 °C (Table 15).

**Table 15. Variation of Extraction Temperature with the AlerTox ELISA Almond Kit.**

Sample	Result 60 °C	Result 25 °C	Result 40 °C	Result 70 °C
Cookie, 0 ppm almond	<b>&lt; 0.4 ppm</b>	< 0.4 ppm	< 0.4 ppm	< 0.4 ppm
Cookie, 4 ppm almond	<b>3.99 ppm</b>	4.17 ppm	4.15 ppm	4.27 ppm

Considering the precision data (Section 2), the results did *not* differ significantly.

## 8.2 Variation of Extraction Time

The extraction time, defined as 15 min, was changed to 10 and 20 min (Table 16).

**Table 16. Variation of Extraction Time with the AlerTox ELISA Almond Kit.**

Sample	Result 15 min	Result 10 min	Result 20 min
Cookie, 0 ppm almond	<b>&lt; 0.4 ppm</b>	< 0.4 ppm	< 0.4 ppm
Cookie, 4 ppm almond	<b>3.72 ppm</b>	4.10 ppm	4.22 ppm

Considering the precision data (Section 2), the results did *not* differ significantly.

## 8.3 Drift

In contrast to the test procedure defined in the instruction manual, the incubation time of the samples was extended and reduced by 4 minutes compared to the calibrators (20 min) (Table 17).

**Table 17. Drift in the AlerTox ELISA Almond Kit.**

Sample	Result 20 min	Result 16 min	Result 24 min
Cookie, 0 ppm almond	<b>&lt; 0.4 ppm</b>	< 0.4 ppm	< 0.4 ppm
Cookie, 4 ppm almond	<b>3.60 ppm</b>	3.07 ppm	4.25 ppm

The results differed significantly. Drift in extensive test runs should be avoided by pipetting calibrators once before and once after the samples and using the mean value for the calculation.

## 9. Summary

- The AlerTox ELISA Almond Kit is sensitive with an LOD of 0.06 ppm (mg of whole almond per kg or L of sample) and an LOQ of 0.5 ppm whole almond (mg/kg or mg/L).
- The AlerTox ELISA Almond Kit provides consistent, precise results, as indicated by results of the intra- and inter-assay studies (mean CV < 10%, Section 2).
- Robust recovery (Section 4) and linearity (Section 5) mean that samples and their dilutions are quantitated accurately.
- Extensive cross-reactivity studies showed high specificity (Section 6) except for apricot kernel (100% cross-reactivity).
- Drift in large test runs can be avoided by pipetting calibrators once before and once after the samples and using the calibrator mean values for the calculations.