

BAX® System MP Media

BAX® System MP Media is a selective enrichment broth for use with the BAX® System PCR assays for detecting both *E. coli* and *Salmonella*. This media allows you to screen samples for both target organisms from a single enrichment. This proprietary medium provides nutrients and stable pH conditions to improve the recovery of *E. coli* and *Salmonella* cells from raw ground beef, beef trim, spinach and lettuce.

Appearance

Bulk Media

Dehydrated: Light-colored, free-flowing powder

Rehydrated: Caramel-colored solution

StatMedia™ Soluble Packets

Dehydrated: Light colored powder in clear film

Rehydrated: Caramel-colored solution, slightly cloudy

Preparation

Bulk Media

Dissolve 22.5 grams of media in 1 liter of distilled water.

Final pH 7.2±0.2 at 25°C. Autoclave at 121°C for 15 minutes.

Do not boil. Broth can be stored for up to 2 weeks at 2-8°C.

StatMedia™ Soluble Packets

Sterilize appropriate volume of distilled or deionized water either by autoclaving or filtering (pore size 0.2 µm) into a sterile container. Warm water to 42°C, then aseptically remove packet from pouch and add to the sterilized water. Mix gently until dissolved. Final pH 7.2±0.2 at 25°C. No need to autoclave if used within 4 hours. Beyond 4 hours, broth should be autoclaved at 121°C for 15 minutes before use. Do not boil. Autoclaved broth can be stored for up to 2 weeks at 2-8°C.



Storage and Stability

Bulk Media

Store dehydrated media in original container, with the lid tightly closed and away from bright light, below 25°C. When stored as directed, bulk media will remain stable until the expiration date stamped on the label of the container.

StatMedia™ Soluble Packets

Store soluble packets in tightly resealed pouches away from moisture at 10-25°C. Use within one year of the manufacturing date stamped on the pouch label.

HYGIENA PRODUCT CODE	DESCRIPTION	QUANTITY
MED2003	BAX® System MP Media	2.5 kg
MED2016	BAX® System StatMedia™ MP Media	20x5x33.75 g
MED2029	BAX® System MP Media	10 kg