

One Health Diagnostics™

Julie Weller, Christine Chapman, Deja Latney, David Luedeke, and Andrew Mason Hygiena<sup>®</sup>, 2 Boulden Circle, New Castle, DE 19720

#### INTRODUCTION:

Prevalence studies have demonstrated a strong association between pre- and post-harvest loads of *Campylobacter* in poultry (1). Regulatory testing of *Campylobacter* in poultry has largely focused on processing plants; however, if an integrated approach was used to analyze both on the farm and processing samples, strategies to control and reduce *Campylobacter* can be improved.

## **PURPOSE:**

The purpose of this study was to evaluate the BAX® System Real-Time PCR assay for the detection of *Campylobacter* species from boot swabs compared to culture.

#### **REGISTERED TRADEMARKS:**

BAX® is a registered trademark of Hygiena for its line of equipment, reagents and software used to analyze samples for microbial contamination.

# Detection of *Campylobacter* from Boot Swabs Using Hygiena's BAX® System Real-Time PCR Assay

BAX<sup>®</sup> System Q7

BAX<sup>®</sup> System X 5

foodproof®

microproof®

### **METHODS:**

Eighteen boot swabs were provided by an industry partner and evaluated for naturally occurring *Campylobacter*. Samples were enriched within 48 hours of collection using 100 mL of prewarmed (37 °C) BPW. After briefly homogenizing, 27.5 mL was transferred to a 50 mL conical tube and 27.5 mL of pre-warmed (42 °C) double-strength Bolton broth with 2X supplement was added. Tubes were tightly capped and incubated aerobically at 42 °C for 48 hours.

Samples were analyzed by real-time PCR and confirmed by culture regardless of presumptive result using modified Charcoal Cefoperazone Deoxycholate Agar (mCCDA) and Campy Cefex Agar.

#### **RESULTS:**

95% CI

Campylobacter was detected by real-time PCR in 9/18 samples at 48 hours. All PCR results matched culture with 100% agreement. See table 1.

Statistical significance using the probability of detection (POD) determined there were no differences between PCR and culture.

Paired Study Results for Boot Swabs	
MPN/Test Portion	0.69
N	18
FN	0
FP	0
Method Agreement	18
∑di	0
dPOD	0
sd	0
SE dPOD	0

-0.15, 0.15


Table 1. BAX System Method

#### Table Legend


MPN/Test Portion = Most Probable Number is based on the POD of BAX System fractional test portions N = Number of test portions FN = Number of false negatives FP = Number of false positives  $\sum d_i = The$  difference of the replicates tested by the BAX method and reference culture method dPOD = Difference between the mean of differences  $s_d = The$  standard deviation of the differences SE dPOD = The standard error of the dPOD does not contain zero, then the difference is statistically significant at the 5% level

# SIGNIFICANCE:

Overall, these results demonstrate the ability of the BAX System to accurately detect *Campylobacter* from boot swabs providing the industry with a way to assess prevalence and implement effective farm management practices.



# FLOW CHART:





# REFERENCES:

1. Berghaus, R. D., Thayer, S. G., Law, B. F., Mild, R. M., Hofacre, C. L., Singer, R. S. 2013. Enumeration of *Salmonella* and *Campylobacter* spp. in Environmental Farm Samples and Processing Plant Carcass Rinses from Commercial Broiler Chicken Flocks. *Appl Environ Microbiol*. 79(13): 4106-4114.