

Technical Bulletin: Detection of *Salmonella* and *Listeria monocytogenes* from Cheese Using the BAX[®] System Real-Time PCR Assays

An unpaired matrix study was conducted at an independent laboratory to measure the performance of the BAX[®] System Real-Time PCR Assay for *Salmonella* and Real-Time PCR Assay for *Listeria monocytogenes* against the United States Food and Drug Administration's reference methods to detect each target organism in semi-hard cheese. Samples tested in this study were artificially inoculated at levels expected to produce low (0.2-2 CFU/test portion) or high (5 CFU/test portion) spike levels after 48-72 hours of equilibration at 4°C. The results obtained for each method were compared using the probability of detection (POD). For both *Salmonella* and *L. monocytogenes*, the BAX[®] System method demonstrated equivalent performance to the reference method.

Introduction

Cheese is a popular ready-to-eat dairy product in many countries. Cheese can be produced in a variety of flavors and textures providing consumers with a rich source of fat, calcium, phosphorus and protein. These nutrients however, also favor the growth of pathogenic bacteria (1). The risk for contamination is the highest for soft cheeses and cheeses made with unpasteurized raw milk. As a result, these high moisture cheeses account for the majority of foodborne outbreaks. Other types of cheese such as semi-hard are also susceptible to contamination and because these types are consumed more often, the microbiological safety of these must also be measured.

Sample Preparation and Enrichment

Salmonella Newport ATCC 6962 and *Listeria monocytogenes* ATCC 19115 were used to inoculate samples of semi-hard cheese independently. For each target organism, cheese was inoculated with a broth culture in bulk to create a low fractional spike level or high spike level. Each master sample was thoroughly mixed to achieve equal distribution of the inoculum and stored at 4°C for 48-72 hours. The master samples were then separated into 25 g test portions to create 20 low-spike and 5 high-spike

samples per method. An additional 5 samples per method were left uninoculated to serve as negative controls.

For the BAX[®] System Real-Time *Salmonella* method, 375 g test portions were analyzed by combining a 25 g sample of inoculated cheese with 350 g of uninoculated cheese. Samples were homogenized with 1500 mL pre-warmed (35°C) Buffered peptone water (BPW) and incubated at 35°C for 22-26 hours. Samples were tested by the BAX[®] System method directly from the primary enrichment and after a BHI regrowth. The regrowth was performed by transferring 10 µL of the primary enrichment to 500 µL pre-warmed (37°C) BHI broth and incubating at 37°C for 3 hours before proceeding to the lysis procedure described below.

For the BAX[®] System Real-Time *Listeria* methods, 375 g test portions were analyzed by combining a 25 g sample of inoculated cheese with 350 g of uninoculated cheese. Samples were homogenized with 1500 mL pre-warmed (20-35°C) 24 LEB Complete media and incubated at 35°C for 26-48 hours. Samples were tested by the BAX[®] System method at 26 and 48 hours.

For the FDA BAM reference methods, 25 g samples were analyzed for both *Salmonella* and *Listeria*. *Salmonella* samples were

homogenized with 225 mL of Lactose broth (LB) and incubated at 35°C for 22-26 hours. *Listeria* samples were homogenized with 225 mL of BLEB containing pyruvate and incubated at 30°C for 4 hours. After 4 hours, solutions of three selective agents (acriflavine, cycloheximide and nalidixic acid) were added to the enrichment, mixed and incubated at 30°C for an additional 44 hours.

Method

BAX® System Method – For *Salmonella* samples, 5 µL of enrichment was added to 200 µL prepared lysis reagent (150 µL of protease to one 12 mL bottle of lysis buffer) in cluster tubes. Lysis was performed by heating tubes for 20 minutes at 37°C and 10 minutes at 95°C, and then cooling tubes at 4°C. Real-Time *Salmonella* PCR tubes were hydrated with 30 µL of lysate and held for 10 minutes on a chilled (2-8°C) cold block. All PCR tubes were then loaded into the BAX® System Q7 instrument, and a full process was run according to the procedure described in the BAX® System User Guide.

For *Listeria* samples, 5 µL of enrichment was added to 200 µL prepared lysis reagent (150 µL of protease and 200 µL of lysing agent 2

to one 12 mL bottle of lysis buffer) in cluster tubes. Lysis was performed by heating tubes for 30 minutes at 55°C and 10 minutes at 95°C, and then cooling tubes at 4°C. Real-Time Genus *Listeria* and Real-Time *L. monocytogenes* PCR tubes were hydrated with 30 µL of lysate. All PCR tubes were loaded into the BAX® System Q7 instrument, and a full process was run according to the procedure described in the BAX® System User Guide.

Reference Method – Each sample was culture confirmed regardless of presumptive BAX® System results following the FDA BAM Chapter 5 for *Salmonella* or the FDA BAM Chapter 10 for *Listeria monocytogenes*.

Results and Discussion

Statistical analysis using the probability of detection (POD) and the difference in POD (dPOD) values were calculated with 95% confidence intervals to compare the results between the BAX® System method presumptive and confirmed results (Table 1) and between the BAX® System method and reference method (Table 2).

Table 1. BAX® System Presumptive vs. Confirmed Results												
Sample Type	Test Method	Target Organism	MPN/25 g	Test Portions	BAX® System Presumptive			BAX® System Confirmed			dPOD _{CP}	95% CI
					X	POD _{CP}	95% CI	X	POD _{CC}	95% CI		
Semi-hard Cheese (375 g)	BAX® System 24 h	<i>Salmonella</i> Newport ATCC 6962	Control	5	0	0.00	0.00, 0.45	0	0.00	0.00, 0.45	0.00	-0.43, 0.43
			1.45	20	13	0.65	0.43, 0.82	14	0.70	0.48, 0.85	-0.05	-0.32, 0.23
			7.28	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.43, 0.43
	BAX® System 26 h	<i>L. monocytogenes</i> ATCC 19115	Control	5	0	0.00	0.00, 0.45	0	0.00	0.00, 0.45	0.00	-0.43, 0.43
			1.02	20	9	0.45	0.26, 0.66	10	0.50	0.30, 0.70	-0.05	-0.33, 0.24
			7.28	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.43, 0.43
	BAX® System 48 h	<i>L. monocytogenes</i> ATCC 19115	Control	5	0	0.00	0.00, 0.45	0	0.00	0.00, 0.45	0.00	-0.43, 0.43
			1.02	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	-0.28, 0.28
			7.28	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.43, 0.43

Results reported for *Salmonella* show data with and without a BHI regrowth.

Results reported for *Listeria monocytogenes* show data from Real-Time Genus *Listeria* and Real-Time *L. monocytogenes*.

MPN/25 g = Most Probable Number is based on the POD of reference method test portions

X = Number of positive test portions

POD_{CP} = BAX® method presumptive positive results divided by the total number of test portions

POD_{CC} = BAX® method confirmed positive results divided by the total number of test portions

dPOD_{CP} = Difference between the BAX® method presumptive result and BAX® method confirmed result POD values

95% CI = If the confidence interval of dPOD does not contain zero, then the difference is statistically significant at the 5%

Table 2. BAX® System Results vs. Reference Method Results

Sample Type	Test Method Comparison	Target Organism	MPN/25 g	Test Portions	BAX® System Method			Reference Method			dPOD _c	95% CI
					X	POD _c	95% CI	X	POD _r	95% CI		
Semi-hard Cheese	BAX® System 24 h / Reference	<i>Salmonella</i> Newport ATCC 6962	Control	5	0	0.00	0.00, 0.45	0	0.00	0.00, 0.45	0.00	-0.43, 0.43
			1.45	20	14	0.70	0.48, 0.85	15	0.75	0.53, 0.88	-0.05	-0.31, 0.22
			7.28	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.43, 0.43
	BAX® System 26 h / Reference	<i>L. monocytogenes</i> ATCC 19115	Control	5	0	0.00	0.00, 0.45	0	0.00	0.00, 0.45	0.00	-0.43, 0.43
			1.02	20	9	0.45	0.26, 0.66	13	0.65	0.43, 0.82	-0.20	-0.45, 0.10
			7.28	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.43, 0.43
	BAX® System 48 h / Reference	<i>L. monocytogenes</i> ATCC 19115	Control	5	0	0.00	0.00, 0.45	0	0.00	0.00, 0.45	0.00	-0.43, 0.43
			1.02	20	10	0.50	0.30, 0.70	13	0.65	0.43, 0.82	-0.15	-0.41, 0.14
			7.28	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.43, 0.43

MPN/25 g = Most Probable Number is based on the POD of reference method test portions

X = Number of positive test portions

POD_c = Confirmed BAX® method positive results divided by the total number of test portions

POD_r = Confirmed reference method positive results divided by the total number of test portions

dPOD_c = Difference between the BAX® method and reference method POD values

95% CI = If the confidence interval of dPOD does not contain zero, then the difference is statistically significant at the 5% level

For *Salmonella* samples enriched in BPW, the BAX® System method returned identical positive results for 13/20 low spiked and 5/5 high spiked samples with and without BHI regrowth. These results matched culture except for 1 additional low spiked sample that screened negative and culture confirmed. The corresponding samples enriched using the FDA BAM method returned culture positive results for 15/20 low spiked and 5/5 high spiked samples. For all method comparisons, no significant statistical difference was observed.

For *Listeria* samples enriched in 24 LEB complete media, the BAX® System method returned positive results for 9/20 low spiked samples at 26 hours and 10/20 low spiked samples at 48 hours. All 5 high spiked samples were positive for both timepoints. All 10 low spiked samples and 5 high spiked samples confirmed as positive. The corresponding samples enriched using the FDA BAM method returned culture positive results for 13/20 low spiked and 5/5 high spiked samples. For all method

comparisons, no significant statistical difference was observed.

Conclusions

These results demonstrate that the BAX® System Real-Time PCR Assays for *Salmonella*, Genus *Listeria* and *L. monocytogenes* performs equivalent to the reference methods for detecting each target organism in semi-hard cheese using the following enrichment protocols:

- For *Salmonella*, homogenize 375 g sample with 1500 mL pre-warmed (35°C) BPW and incubate at 35°C for 24 hours.
- For *Listeria*, homogenize 375 g sample with 1500 mL pre-warmed (20-35°C) 24 LEB Complete media and incubate at 35°C for 26-48 hours.

References

1. Kousta, M., Mataragas, M., Skandamis, P., Drosinos E. H. 2010. Prevalence and sources of cheese contamination with pathogens at farm and processing levels. Food Control. 21 (6): 805-815.