


Hygiena's Positive Control Kit

The Positive Control Kit validates the efficacy and stability of Hygiena® ATP test devices. The kit comes with 25 sealed glass vials containing lyophilized adenosine triphosphate (ATP) and sugars to provide predictable results when sampling with test devices.

Benefits:

- Validates performance of Hygiena test devices
- No lab equipment needed
- No reconstitutions or complicated sample preparation
- Helps to meet GMP and other quality program requirements



Compatible Products

- UltraSnap® ATP Surface Test
- SuperSnap® High Sensitivity ATP Surface Test
- SnapShot Universal ATP Surface Test
- SpotCheck Plus Food Safety Residue Test

Catalog No.	Description	Quantity
CK25	ATP Positive Control Kit	1
US2020	UltraSnap Surface ATP Test	100
SUS3000	SuperSnap High Sensitivity ATP Surface Test	100
See catalog	SnapShot Universal ATP Test	100
GL100	SpotCheck Plus	100

See an instructional demo at [YouTube HygienaTV](#) or [Twitter HygienaTV](#)