

Commodities validated for Hygiena™ Mycotoxin ELISA Test Kits • May 2020



Commodity	Aflatoxin B1	Aflatoxin M1 Sensitive	T-2 Toxin	Aflatoxin M1 Urine	Aflatoxin M1 Ultra	Total Aflatoxin	Total Aflatoxin LM	Mycotox Total Aflatoxin	Deoxynivalenol	Fumonisin	Fumonisin Urine	Aflatoxin B1 Low Matrix	Ochratoxin A Universal	Zearalenone Low Matrix
Almond							X					X		
Animal Feed			X				X		X			X		X
Barbecue Sauce							X					X		
Barley									X					
Beer													X	
Cannabis Flower							X					X	X	
Cardamom							X					X		
Cereals													X	X
Char Siu Sauce							X					X		
Chickpeas							X					X		
Chilli Powder							X					X		
Chilli Sauce							X					X		
Cocoa													X	
Cocoa Butter													X	
Coconut							X					X		
Coffee													X	
Coriander							X					X		
Corn	X		X			X	X	X	X	X		X	X	
Corn Oil							X					X		
Cottonseed	X					X	X					X		
Curry Sauce							X					X		
Dumpling Sauce							X					X		
Flaxseed							X					X		
Ginger							X					X		
Grape Juice													X	
Grape Must													X	
Hay							X					X		
Hemp Seed							X					X	X	
Hoisin Sauce							X					X		
Hunan Sauce							X					X		
Infant Formula							X					X		

Commodity	Aflatoxin B1	Aflatoxin M1 Sensitive	T-2 Toxin	Aflatoxin M1 Urine	Aflatoxin M1 Ultra	Total Aflatoxin	Total Aflatoxin LM	Mycotox Total Aflatoxin	Deoxynivalenol	Fumonisin	Fumonisin Urine	Aflatoxin B1 Low Matrix	Ochratoxin A Universal	Zearalenone Low Matrix
Kolo (grain)							X					X		
Macadamia							X					X		
Milk													X	
Milk Powder, Skim		X			X									
Milk, Raw		X			X									
Milk, Whole Hmgnzd		X			X									
Nutmeg							X					X		
Oyster Sauce							X					X		
Paprika							X					X		
Peanut Oil							X					X		
Peanut							X					X		
Pecans							X					X		
Pistachio							X					X		
Plum Sauce							X					X		
Rice							X					X		
Safflower Oil							X					X		
Serum													X	
Sesame Oil							X					X		
Sesame Seed							X					X		
Shiro (spice)							X					X		
Shrimp Paste							X					X		
Snaplage							X					X		
Soy Sauce							X					X		
Soybean							X					X		
Spices													X	
Sunflower Seed							X					X		
Urine				X							X			
Vegetable Oil							X					X		
Wheat			X				X		X			X		
Wheat Bran									X					
Wine													X	
Yogurt					X									