

Hygiena™ Biomass Kit


The Hygiena™ Biomass Kit is a high-precision test for measuring total and free ATP in liquid samples, offering expanded flexibility over all-in-one AquaSnap™ Water ATP test devices.

Benefits:

- Detects microbial and non-microbial contamination
- Uses patented liquid-stable reagent for better accuracy and reproducibility
- Enables cleanup with high-quenching particulates or fibers in aqueous samples
- 12 month shelf-life at refrigerated temperature (2 - 8 °C)
- ATP standard for reagent activity testing, and 100 cuvettes included.

Additional Requirements:

- Luminometer (EnSURE™ or SystemSURE Plus)
- Hold-Rite applicator (HR002)
- Pipette and ATP free tips for 100 µl sample volume



Catalog No.	Description	Quantity
CCK-2	Biomass Kit for Celsis™ instruments	1 kit
CCK-4	Biomass Kit for Hygiena SystemSURE II, SystemSURE Plus, and EnSURE™ luminometers	1 kit
HR002	Hold-Rite Applicator	1



Find support documents, instructional videos, and more at www.hygiena.com