



Transport Swab


Environmental Sample
Collection and Dilution Device

Hygiena's Transport Swab is the ideal sampling device for customers looking for an environmental sample collection device containing more than 1 mL of enrichment media. Transport Swab makes sample collection efficient and easy.

Benefits:

- Ready-to-use and convenient—pre-filled bottles come sterilized and have ample head space for additional sample
- All-in-one sample collection system
- More broth for high-volume samples
- Available in 4 mL and 10 mL formats
- Available with Lethen or Buffered Peptone Water
- Perfect for plates or dried media film
- Collect your samples by swabbing the surface and recapping

Transport Swab Procedure



Catalog No.	Description	Quantity
HTS-4	Transport Swab (Lethen 4 mL)	100
HTS-10	Transport Swab (Lethen 10 mL)	100
HTS-10BPW	Transport Swab (Buffered Peptone Water 10 mL)	100