

	KIT5001	KIT5003	KIT5004	KIT5005	KIT5006	KIT5007	KIT5008	KIT5009	KIT5010	KIT5014	KIT5015	KIT5017	KIT5018
Commodity	Aflatoxin M1	Aflatoxin M1 ULTRA	Aflatoxin B1 Rapid	Aflatoxin B1 Low Matrix	Total Aflatoxin Low Matrix	Total Aflatoxin Rapid	MycoTox Total Aflatoxin	Aflatoxin Total Hydro	Deoxynivalenol	Fumonisin Hydro	Ochratoxin A Universal	T-2 Toxin	Zearalenone Low Matrix
Kolo (grain)				X	X								
Macadamia				X	X								
Milk											X		
Milk Powder: Skim	X												
Milk: Raw	X												
Milk: Whole Hmgznd	X												
Nutmeg				X	X								
Oyster Sauce				X	X								
Paprika				X	X								
Peanut Oil				X	X								
Peanut				X	X								
Pecans				X	X								
Pistachio				X	X								
Plum Sauce				X	X								
Rice				X	X								
Safflower Oil				X	X								
Serum											X		
Sesame Oil				X	X								
Sesame Seed				X	X								
Shiro (spice)				X	X								
Shrimp Paste				X	X								
Snaplage				X	X								
Soy Sauce				X	X								
Soybean				X	X								
Spices											X		
Sunflower Seed				X	X								
Vegetable Oil				X	X								
Wheat				X	X				X			X	
Wheat Bran									X				
Wine											X		
Yogurt		X											