

Stick Sponge™ | Sample Collection




Hygiena™ Stick Sponges are biocide-free cellulose sponges attached to a plastic handle and supplied in a sterile, tie-off bag. The easy-to-release stick makes it easier to reach into drains, pipes, and around equipment.

Benefits:

- For surface sampling and hard to reach areas
- 1.5" x 3" biocide-free cellulose sponge
- Sterilized by gamma irradiation
- Pre-moistened with a variety of diluents
- Sample collection without directly touching the sponge
- Easy-to-open tie-off bag
- Large write-on area for easy sample identification
- 12 - 18-month shelf life for hydrated sponges
- Attached sterile gloves optional

Stick Sponge Procedure



See an instructional demo at [HygienaTV](#) or [HygienaTV](#)

Catalog No.	Description	Quantity
SS100BPW	Stick Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 10 bags of 10	100
SS100BPWG	Stick Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 10 bags of 10 with gloves	100
SS100NB	Stick Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 10 bags of 10	100
SS100NBG	Stick Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 10 bags of 10 with gloves	100
SS100LET	Stick Sponge, 1.5" x 3" Sponge w/10mL Lethen Broth, 10 bags of 10	100
SS100LETG	Stick Sponge, 1.5" x 3" Sponge w/10mL Lethen Broth, 10 bags of 10 with gloves	100



Find support documents, instructional videos, and more at www.hygiena.com

Sponge'n Bag™ | Sample Collection




Hygiena™ Sponge'n Bags are biocide-free cellulose sponges supplied in a sterile, tie-off bag. The sponge provides a simple and convenient device for collecting environmental samples from large surface areas. Sponges may be ordered pre-moistened with a variety of different diluents.

Benefits:

- 1.5" x 3" biocide-free cellulose sponge
- Sterilized by gamma irradiation
- Pre-moistened with a variety of diluents
- Easy-to-open tie-off bag
- Large write-on area for easy sample identification
- 12 - 18-month shelf life for hydrated sponges
- Attached sterile gloves optional

Sponge'n Bag Procedure



See an instructional demo at [HygienaTV](#) or [HygienaTV](#)

Catalog No.	Description	Quantity
HSB100BPW	Sponge'n Bag, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 10 bags of 10	100
HSB100BPWG	Sponge'n Bag, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 10 bags of 10 with gloves	100
HSB100NB	Sponge'n Bag, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 10 bags of 10	100
HSB100NBG	Sponge'n Bag, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 10 bags of 10 with gloves	100
HSB100LET	Sponge'n Bag, 1.5" x 3" Sponge w/10mL Lethen Broth, 10 bags of 10	100
HSB100LETG	Sponge'n Bag, 1.5" x 3" Sponge w/10mL Lethen Broth, 10 bags of 10 with gloves	100
HSB100DENB	Sponge'n Bag, 1.5" x 3" Sponge w/10mL D/E Neutralizing Broth, 10 bags of 10	100
HSB100DENBG	Sponge'n Bag, 1.5" x 3" Sponge w/10mL D/E Neutralizing Broth, 10 bags of 10 with gloves	100



Find support documents, instructional videos, and more at www.hygiena.com