

**Safety data sheet**  
**according to Regulation (EC) No 1907/2006, Article 31**

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

---

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****· 1.1 Product identifier****· Trade name: foodproof GMO Kits, ready-to-use****· Product variants:**

KIT230016 / foodproof® CaMV Detection Kit (Ref. F 302 39)

KIT230019 / foodproof® GMO A2704-12 Soya Quantification Kit (Ref. F 302 66)

KIT230053 / foodproof® GMO RR 2 Yield Soya Quantification Kit (Ref. R 302 35)

**· 1.2 Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

**· Application of the substance / the mixture**

PCR-Kit for the detection of certain DNA. Target DNA as specified on separate enclosure.

**· 1.3 Details of the supplier of the safety data sheet****· Manufacturer/Supplier:**

Hygiena Diagnostics GmbH

Hermannswerder 17

14473 Potsdam

GERMANY

Phone: +49 (0) 331-23 00 200

www.hygiena.com

BIOTECON Diagnostics GmbH

Hermannswerder 17

14473 Potsdam

GERMANY

Phone: +49 (0) 331-23 00 200

www.bc-diagnostics.com

**· Informing department:**

Phone: +49 (0) 331-23 00 200

Fax: +49 (0) 331-23 00 299

**· 1.4 Emergency telephone number:**

Phone: +49 (0) 331-23 00 200

Office hours: 09:30 - 16:00

---

**\* SECTION 2: Hazards identification****· 2.1 Classification of the substance or mixture****· Classification according to Regulation (EC) No 1272/2008**

The product is not classified, according to the CLP regulation.

**· Additional information:** Set composed of separate, different vials.**· 2.2 Label elements****· Labelling according to Regulation (EC) No 1272/2008** Void**· Hazard pictograms** Void**· Signal word** Void**· Hazard statements** Void**· Additional information:**

Please refer to further labelling elements in section 15 of this safety data sheet, if applicable.

**· 2.3 Other hazards****· Results of PBT and vPvB assessment****· PBT:** Not applicable.

(Contd. on page 2)

# Safety data sheet

## according to Regulation (EC) No 1907/2006, Article 31

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

**Trade name: foodproof GMO Kits, ready-to-use**

(Contd. from page 1)

- **vPvB:** Not applicable.
- **Determination of endocrine-disrupting properties**  
The product does not contain substances with endocrine disrupting properties.

### SECTION 3: Composition/information on ingredients

- **3.2 Mixtures**
  - **Description:** PCR-Kit for the detection of certain DNA. Target DNA as specified on separate enclosure.
  - **Dangerous components:** Void
  - **Additional information** For the wording of the listed hazard phrases refer to section 16.

### SECTION 4: First aid measures

- **4.1 Description of first aid measures**
  - **General information**  
Remove contaminated clothing immediately.  
No special measures required.
  - **After inhalation** Supply fresh air; consult doctor in case of symptoms.
  - **After skin contact**  
Wash with water and soap.  
If skin irritation continues, consult a doctor.
  - **After eye contact**  
Rinse opened eye for several minutes under running water.  
Remove contact lenses if possible.  
In case of permanent aches and pains please go and see the doctor.
  - **After swallowing** Rinse out mouth and then drink plenty of water.
- **4.2 Most important symptoms and effects, both acute and delayed**  
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
  - **Suitable extinguishing agents**  
The product itself does not burn.  
Use fire fighting measures that suit the environment.
  - **For safety reasons unsuitable extinguishing agents** None known.
- **5.2 Special hazards arising from the substance or mixture**  
Formation of toxic gases is possible during heating or in case of fire.
- **5.3 Advice for firefighters**
  - **Protective equipment:** Wear self-contained breathing apparatus.
  - **Additional information**  
Collect contaminated fire fighting water separately. It must not enter drains. Provide sufficient fire fighting water retention.  
Dispose of fire debris and contaminated fire fighting water in accordance with official regulations.

### SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**  
Avoid contact with the product.
- **6.2 Environmental precautions:** Do not allow to enter drainage system, surface or ground water.
- **6.3 Methods and material for containment and cleaning up:**  
Send for recovery or disposal in suitable containers.

(Contd. on page 3)

# Safety data sheet

## according to Regulation (EC) No 1907/2006, Article 31

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

---

**Trade name: foodproof GMO Kits, ready-to-use**


---

(Contd. from page 2)

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders).  
Dispose of the material collected according to regulations.

- **6.4 Reference to other sections**

- See Section 7 for information on safe handling

- See Section 8 for information on personal protection equipment.

- See Section 13 for information on disposal.

---

## SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**

- Store in cool, dry place in tightly closed containers.

- Keep away from heat and direct sunlight.

- Avoid contact with eyes and skin.

- No special measures required.

- **Information about protection against explosions and fires:** No special measures required.

- **Handling**

- Do not eat, drink or smoke while working.

- Wash hands during breaks and at the end of the work.

- **7.2 Conditions for safe storage, including any incompatibilities**

- **Storage**

- **Requirements to be met by storerooms and containers:** Store only in the original container.

- **Information about storage in one common storage facility:**

- Keep away from strong oxidizing, alkalis and acidic materials.

- **Further information about storage conditions:**

- Store in the dark.

- Protect from the effects of light.

- **Recommended storage temperature:** < 25°C.

- **7.3 Specific end use(s)** No further relevant information available.

---

## SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**

- **Components with limit values that require monitoring at the workplace:**

- The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- **Additional information:** The lists that were valid during the compilation were used as basis.

- **8.2 Exposure controls**

- **Individual protection measures, such as personal protective equipment**

- **Breathing equipment:** Not required when used as intended.

- **Hand protection**

- Protective gloves are recommended at frequent and/or long-term skin contact with the product.

- To avoid skin problems reduce the wearing of gloves to the required minimum.

- **Material of gloves**

- The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

- The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **For the permanent contact in work areas without heightened risk of injury (e.g. Laboratory) gloves made of the following material are suitable:**

- Nitrile rubber, NBR

- **Eye/face protection** Safety glasses recommended if possibility of liquid splashes exists.

(Contd. on page 4)

# Safety data sheet

## according to Regulation (EC) No 1907/2006, Article 31

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

**Trade name: foodproof GMO Kits, ready-to-use**

(Contd. from page 3)

- **Body protection:** Wear suitable protective clothing.

## SECTION 9: Physical and chemical properties

### · 9.1 Information on basic physical and chemical properties

#### · General Information

- |   |                              |
|---|------------------------------|
| · <b>Physical state</b>   | Liquid                       |
| · <b>Colour:</b>  | Yellowish                    |
| · <b>Odour:</b>   | Light                        |
| · <b>Odour threshold:</b>   | Not determined.              |
| · <b>Melting point/freezing point:</b>                            | Not determined               |
| · <b>Boiling point or initial boiling point and boiling range</b> | Not determined               |
| · <b>Flammability</b>   | Not applicable.              |
| · <b>Lower and upper explosion limit</b>                          |                              |
| · <b>Lower:</b>   | Not determined.              |
| · <b>Upper:</b>   | Not determined.              |
| · <b>Flash point:</b>   | Not applicable               |
| · <b>Auto-ignition temperature:</b>                               | Product is not selfigniting. |
| · <b>Decomposition temperature:</b>                               | Not determined.              |
| · <b>SADT</b>   |                              |
| · <b>pH</b>   | Not determined.              |
| · <b>Viscosity:</b>   |                              |
| · <b>Kinematic viscosity</b>                                      | Not determined.              |
| · <b>dynamic:</b>   | Not determined.              |
| · <b>Solubility</b>   |                              |
| · <b>Water:</b>   | mixable                      |
| · <b>Partition coefficient n-octanol/water (log value)</b>        | Not determined.              |
| · <b>Vapour pressure:</b>   | Not determined.              |
| · <b>Density and/or relative density</b>                          |                              |
| · <b>Density</b>  | Not determined               |
| · <b>Relative density</b>   | Not determined.              |
| · <b>Vapour density</b>   | Not determined.              |

### · 9.2 Other information

- |  |                              |
|--|------------------------------|
| · <b>Appearance:</b>   |                              |
| · <b>Form:</b>   | Fluid                        |
| · <b>Important information on protection of health and environment, and on safety.</b> |                              |
| · <b>Self-inflammability:</b>  | Product is not selfigniting. |
| · <b>Explosive properties:</b>   | Product is not explosive.    |
| · <b>Change in condition</b>   |                              |
| · <b>Evaporation rate</b>  | Not determined.              |

### · Information with regard to physical hazard classes

- |  |      |
|--|------|
| · <b>Explosives</b>                            | Void |
| · <b>Flammable gases</b>                       | Void |
| · <b>Aerosols</b>                              | Void |
| · <b>Oxidising gases</b>                       | Void |
| · <b>Gases under pressure</b>                  | Void |
| · <b>Flammable liquids</b>                     | Void |
| · <b>Flammable solids</b>                      | Void |
| · <b>Self-reactive substances and mixtures</b> | Void |
| · <b>Pyrophoric liquids</b>                    | Void |
| · <b>Pyrophoric solids</b>                     | Void |

(Contd. on page 5)

# Safety data sheet

## according to Regulation (EC) No 1907/2006, Article 31

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

Trade name: foodproof GMO Kits, ready-to-use

(Contd. from page 4)

· <b>Self-heating substances and mixtures</b>	Void
· <b>Substances and mixtures, which emit flammable gases in contact with water</b>	Void
· <b>Oxidising liquids</b>	Void
· <b>Oxidising solids</b>	Void
· <b>Organic peroxides</b>	Void
· <b>Corrosive to metals</b>	Void
· <b>Desensitised explosives</b>	Void

## SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
  - **Thermal decomposition / conditions to be avoided:**  
No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:**  
None in case of intended use and storage in compliance with instructions.

## SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
  - **Acute toxicity** Based on available data, the classification criteria are not met.
  - **Primary irritant effect:**
    - **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
    - **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
  - **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**

### · **Endocrine disrupting properties**

None of the ingredients is listed.

## SECTION 12: Ecological information

- **12.1 Toxicity**
  - **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** No further relevant information available.
- **Other information:** There are no data available about the preparation.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
  - **PBT:** Not applicable.
  - **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties**  
The product does not contain substances with endocrine disrupting properties.
- **12.7 Other adverse effects** No further relevant information available.

(Contd. on page 6)

# Safety data sheet

## according to Regulation (EC) No 1907/2006, Article 31

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

**Trade name: foodproof GMO Kits, ready-to-use**

(Contd. from page 5)

- **Additional ecological information:**

- **General notes:**

- Do not allow undiluted product or large quantities of it to reach ground water, water bodies or sewage system.

## SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**

- **Recommendation**

- Smaller quantities can be disposed with household waste.

- The waste code numbers mentioned are recommendations based on the probable use of the product.

- **Uncleaned packagings:**

- **Recommendation:**

- Dispose of packaging according to regulations on the disposal of packagings.

- Non contaminated packagings can be used for recycling.

- Packagings that cannot be cleaned are to be disposed of in the same manner as the product.

## SECTION 14: Transport information

<ul style="list-style-type: none"> <li>· <b>14.1 UN number or ID number</b></li> <li>· <b>ADR/RID, IMDG, IATA</b></li> </ul>	Void
<ul style="list-style-type: none"> <li>· <b>14.2 UN proper shipping name</b></li> <li>· <b>ADR/RID, IMDG, IATA</b></li> </ul>	Void
<ul style="list-style-type: none"> <li>· <b>14.3 Transport hazard class(es)</b></li> <li>· <b>ADR/RID, ADN, IMDG</b></li> <li>· <b>Class</b></li> </ul>	Void
<ul style="list-style-type: none"> <li>· <b>IATA</b></li> <li>· <b>Class</b></li> </ul>	Void Not Restricted.
<ul style="list-style-type: none"> <li>· <b>14.4 Packing group</b></li> <li>· <b>ADR/RID, IMDG, IATA</b></li> </ul>	Void
<ul style="list-style-type: none"> <li>· <b>14.5 Environmental hazards:</b></li> <li>· <b>Marine pollutant:</b></li> </ul>	No
<ul style="list-style-type: none"> <li>· <b>14.6 Special precautions for user</b></li> </ul>	Not applicable.
<ul style="list-style-type: none"> <li>· <b>14.7 Maritime transport in bulk according to IMO instruments</b></li> </ul>	Not applicable.
<ul style="list-style-type: none"> <li>· <b>Transport/Additional information:</b></li> </ul>	Not dangerous according to the above specifications.
<ul style="list-style-type: none"> <li>· <b>UN "Model Regulation":</b></li> </ul>	Void

## SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **Directive 2012/18/EU**

- **Named dangerous substances - ANNEX I** None of the ingredients is listed.

(Contd. on page 7)



# Safety data sheet

## according to Regulation (EC) No 1907/2006, Article 31

Printing date 10.02.2025

Version number 2 (replaces version 1)

Revision: 10.02.2025

**Trade name: foodproof GMO Kits, ready-to-use**

(Contd. from page 6)

· **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II**

None of the ingredients is listed.

· **REGULATION (EU) 2019/1148**

· **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))**

None of the ingredients is listed.

· **Annex II - REPORTABLE EXPLOSIVES PRECURSORS**

None of the ingredients is listed.

· **Regulation (EC) No 273/2004 on drug precursors**

None of the ingredients is listed.

· **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors**

None of the ingredients is listed.

· **National regulations**

· **Water hazard class:**

· **Substances of very high concern (SVHC) according to REACH, Article 57**

None of the ingredients is contained.

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

This data sheet considers labelling requirements of Regulation (EU) 2022/692 adapting Regulation (EC) No 1272/2008 (18th ATP of the CLP Regulation).

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

This safety data sheet meets the requirements of Regulation (EU) 2015/830 and 2020/878 amending Annex II of Regulation (EC) 1907/2006.

· **Relevant phrases**

The phrases specified here are no labelling elements for the product but repeat the properties of the ingredients from section 3.

· **Department issuing data specification sheet:**



This Safety Data Sheet has been drawn up in cooperation with:  
DEKRA Assurance Services GmbH, Hanomagstr. 12, D-30449 Hanover, Germany,  
phone: (+49) 511 42079 - 0, reach@dekra.com.

© DEKRA Assurance Services GmbH. Changing this documents is subject to explicit acceptance by DEKRA Assurance Services GmbH.

· **Date of previous version:** 08.10.2023

· **Version number of previous version:** 1

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonised System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service (division of the American Chemical Society)  
PBT: Persistent, Bioaccumulative and Toxic  
SVHC: Substances of Very High Concern  
vPvB: very Persistent and very Bioaccumulative

· **\* Data compared to the previous version altered.**